

# Portable Tactile & Auditory

## TACTILE AND AUDITORY EMDR IN ONE SMALL PACKAGE



*Deluxe*  
**Tac/AudioScan™**



*Advanced*  
**Tac/AudioScan™**



*Basic*  
**Tac/AudioScan™**

Now, you have 3 portable auditory and tactile EMDR instruments to choose from, the **Deluxe Tac/AudioScan**, the **Advanced Tac/AudioScan** and the original **Basic Tac/AudioScan**. All 3 products provide easy to use tactile and

auditory functions for effective EMDR administration. Auditory and tactile modes can be used independently or use both functions together for faster and deeper processing.

**Auditory stimulus:** All models provide alternating tones to each ear via light-weight stereo headphones. The *Deluxe and Advanced Tac/AudioScans* also give you the choice of 3 additional sounds; alternating click, alternating double click (heart beat), and alternating arcade sound (kids love this one). Plus, the *Deluxe & Advanced Tac/AudioScans* also have the capability of connecting to any music source, such as a tape, CD, iPod or MP3 player as shown to the right, turning **your** music or **any** recorded sound into bilateral sound. The best part is that the auditory and tactile functions are synchronized.

**Tactile stimulus:** All models feature fully adjustable tactile stimulation. Tactile intensity is infinitely adjustable from completely off to full maximum strength. The included oval shaped tactile pulsers are extremely small, only 1.5" x 1" x 1/2". They are easy to hold in the hands or position in other convenient places such as in pockets or socks. Additionally, the pulsers may be attached to wrists, arms or legs using commonly available velcro straps or exercise bands (not included).

**Digital display:** The Deluxe model also includes a digital display for count/setup information and a low battery indicator.

**Easy to use:** All controls, adjustments and visual speed indicators are on the front of the unit allowing easy access while the headphone and tactile pulsers connect on the top. Since the auditory and tactile stimulations are synchronized, the speed adjustment regulates both tactile and audio rates. Volume, speed and tactile intensity adjustments have calibrated dials allowing the clinician to record each clients' preferred settings.

**Order yours today:** All models come complete with headphones, standard tactile pulsers, carrying case and battery. The *Deluxe & Advanced* models also include free AC adapter and audio cable. An optional AC adapter is available for Basic Tac/AudioScan. An optional CD player is available for the *Deluxe & Advanced Tac/AudioScan*. With hands-free administration of EMDR, the clinician is better able to closely observe the client's behavior and take notes. All Tac/AudioScan™ models are effective therapeutic tools and are only available from NeuroTek.

™ Trademark of NeuroTek Corporation.

©Copyright 2014 NeuroTek Corporation All rights reserved

Neurotek products are only sold to therapists who have completed an EMDRIA approved training program.



**Tac/AudioScan Feature Table**

	Tactile	Number of	External	Digital	Low Bat.	AC Adapter
Model	Pulsers	Sounds	Auditory	Display	Indicator	Included
Basic	Yes	1	No	No	No	No
Advanced	Yes	4	Yes	No	No	Yes
Deluxe	Yes	4	Yes	Yes	Yes	Yes

**www.neurotekcorp.com**